Delta4 – A New IMRT QA Device

Görgen Nilsson - ScandiDos AB

Purpose

The basic capabilities of a 4D dosimetry and quality control device were evaluated in a clinical setting. QA measurements of dynamic treatment modafties such as IMRT, 4DRT, ARC, IMAT, gating and TomoTherapy were perform. Dose delivered was measured in short time intervals allowing for a comparison of the total dose delivered with the treatment plan

Method and Materials

Simultaneous measurements were performed with a Detai* Oft-device of absolute dose and lose distribution in 4-dimensions including time. Companion was done to a Some of the basic characteristics was done in flar phartom. Some of the basic characteristics was done in flar phartom. Monstroining units and segment weighting in MRT, felsts, beams and composite plan worlfield.



Measurement accuracy using Delta4

Deltat uses a new p-Si diodes. Previous Si-diodes on the market are known to have very good spatial resolution and simple to use but the drawback of over response to low energy photons. The clinical relevance is possible discrepancy in the response at various depths, field sizes and outside the primary field. The new p-Si diode was compared to a Farmer type chamber e.g. NE 257 is a clinical relevant conditions.

DPPD test: Constant build-up, back scatter and field size at detector; IC corrected for recombination, various SDD was used to change dose per pulse without changing energy spectrum.

Field dependency test: Constant SDD; changing field size at various depths; normalization at 10x10 cm2


Depth dependency test: Constant SDD and field size 10x10 cm2; various depths and back scatter total thickness 22cm (same as Delta"), normalization at reference depth



	Energy	Depth (mm)	50x50 m
	6MV	15	99.7
0.00	(6.0)	45	99.7
13		105	99.1
		205	99.1
	10		
	15MV	15	99.8
		45	99.7
	1	105	99.7
- The state of the	n.	205	99.7
		Depth d	
	10v10 em	2	Datie F


Field size dependency Delta ⁴					
Energy		Ratio Delta ⁴ / IC chamber (%)			
	Depth (mm)	50x50 mm ²	100x100	150x150	200x200
6MV	15	99.7	100.0	100.7	101.2
	45	99.7	100.0	101.2	102.2
	105	99.1	100.0	101.7	102.8
	205	99.1	100.0	101.4	102.5
15MV	15	99.8	100.0	99.7	99.9
	45	99.7	100.0	100.4	101.1
	105	99.7	100.0	100.6	101.7
	205	99.7	100.0	100.7	101.3

Depth dependency Delta⁴					
10x10 cm2	Ratio	Ratio Delta ⁴ / IC chamber (%)			
Energy	15 mm	45 mm	105 mm	205 mm	
6MV	99.8	100.0	100.9	100.4	
15MV	101.4	100.0	100.5	100.1	



Out of field dependency Delta ⁴				
		Ratio Delta / IC chamber (%)		
Energy	Distance to field edge (mm)	50x50 mm ²	100x100 mm ²	150x150 mm²
6MV'	20	0.5	0.9	1.3
	40	0.3	0.8	1.0
	60		0.6	0.8
	80			0.7
15MV	20	0.3	0.2	1.1
_	40	0.2	0.3	0.8
	60		0.2	0.7
	80			0.8







QA in various applications using Delta4

Date I measure the done distribution is a modern ear, synchronized with the societative trip place, all descrips panels. Date is strokedually for each follow place together instruction external systems or, add. Existing, earlying, earling, 408T etc. making the system externey's facilities to various treatment techniques. Date from the same measurement was compared to plan date both per fraction, bears and segment (MLC compared note to non available date from TPs) in the too measuring places.

When the system externey's facilities are successful to the system and compared with TPs dates. DVH where compared for the semi-measured dates with the TPS date for the patient structures applied to the phartom to evaluate the significance of and provision.

Measurement and TPS data in the two detector planes displayed as dose deviation (could be DTA, gamma index etc as well). A profile true the 2D array displays 1D data comparing measured data and TPS data.





The State of the Debt Out-driven showed a very small field and depth dependency of less than 1.5%, a decrease in sensitivity of less. The p-St doubte in the Debt Out-driven showed a very small field and depth dependency of less than 1.5%, a decrease in sensitivity of less and the distinguish deviations in critical organs and tumor from other less sensitive tissue.

Conclusion

The Dated measures each often guide individually making in possible to view date in any time-increment as plan, beam, segment, correct good, respiratory shape, and angle increment. Measurements were directly compared with the does destibution of a composite beatment, sowing extensive time in evaluation of good plans. Dates analysis can be referred to beam and sub-hoars level using the dates acquired during the Chi-forcite plane of the country of the co